

Barrel Shaker VKS 75 F control



Description

Big heavy-duty shaker for high loads also during continuous operation

(floor installation only!)

Basic Equipment

Basic device **without shaking plate, without rack system**

With digital display and programmable control

Loading Capacity

Horizontal or vertical barrels up to 80 l capacity

Option:

[Control interface](#)

[USB interface](#)

Details

Motion: reciprocating (horizontal)

Shaking platform: 760 x 600 mm

Max. load: 75 kg

Shaking speed: up to 80 rpm, steplessly adjustable

Stroke: 120 -150 mm in 10 mm steps
Timer: programmable / continuous
Order no.: 6212 000
Accessory: Shaking plate VKS (1-storey)
Rack system Barrel Shaker

Technical specifications

Electrical supply: 100 - 240 V, 50 / 60 Hz please indicate in case of order
Enclosure protection: IP 21
Heat emission: max. 45 W
Ambient temperature: 5°C to 50°C
Relative humidity: ~ 85 %
Dimensions (w x d x h): 1050 x 835 x 250 mm
Weight: 110 kg
Total weight (Shaker + shaking plate + Rack system): 132 kg