

Compact Shaker KS 15 A control



Description

For reproducible shaking tasks

Basic Equipment

Basic device without shaking plate, **without rack system**
With digital display and programmable control.

Loading Capacity

Erlenmeyer flasks	100 ml	20 pieces
Erlenmeyer flasks	250 ml	12 pieces
Erlenmeyer flasks	500 ml	6 pieces
Erlenmeyer flasks	1000 ml	4 pieces
Separating funnels	250 ml	3 pieces

Option:

[Control interface](#)

[USB interface](#)

Details

Motion: orbital

Shaking platform: 400 x 300 mm

Max. load: 15 kg
Shaking speed: 30 - 420 rpm, steplessly adjustable
Stroke: 17 mm
Runtime: programmable / continuous
Order No.: 6170 000
Accessory: Rack System Combifix KS
Universal tray KS
Incubator Hood TH 15

Technical specifications

Electrical supply: 230 V or 115 V, 50 / 60 H please indicate in case of order
Enclosure protection: IP 21
Heat emission: approx. 10 - 15 W
Ambient temperature: 5°C to 50°C
Relative humidity: ~ 85 %
Dimensions (w x d x h): 510 x 490 x 150 mm
Weight: 19 kg