

Hot Plate Shaker SM 30 AT control



Description

Table-top shaker with heating plate and programmable shaking tasks for mixing of difficult media

Basic equipment

Basic device with heating plate, **without rack system**
With digital display and programmable control

Loading Capacity

Erlenmeyer flasks	100 ml	42 pieces
Erlenmeyer flasks	250 / 300 ml	20 pieces
Separating funnels	250 ml	6 pieces
Separating funnels	1000 ml	4 pieces

Option:

[Control interface](#)

[USB interface](#)

Details

Motion:	orbital
Shaking platform:	560 x 400 mm
Max. load:	25 kg

Shaking speed: 15 - 300 rpm, steplessly adjustable
Stroke: 26 mm
Runtime: programmable / continuous
Temperature range: up to 120°C
Control accuracy: approx. +/- 1°C
Order no.: 6200 000
Accessory: Rack System Combifix SM A
Rack System Combifix SM B
Rack System Combifix SM C

Technical specifications

Electrical supply: 230 V, 50 / 60 Hz 115 V version on request
Enclosure protection: IP 21
Max. heating power: 400 W
Ambient temperature: 5°C to 50°C
Relative humidity: ~ 85 %
Dimensions (w x d x h): 680 x 610 x 160 mm
Weight: 40 kg