

Incubator Hood TH 15



Description

Constant temperatures; also for high vessels (e.g. 1000 ml Erlenmeyer flasks).
Can be combined with all compact shakers or used as independent unit

Basic Equipment

Coated metal housing with 3 windows, the front door opens upwards

Safety cut-off of the shaker if combined with a Bühler Compact Shaker
Standard temperature controller 5°C above ambient to +50°C

Option

[Extended temperature range up to +60°C](#)

Details

Temperature range: 5°C above ambient to +50°C
Temperature accuracy: + / - 1% of the set value

Inner dimensions (w x d x h): 90 x 410 x 310 mm (height sufficient for 1000 ml Erlenmeyer flasks)

Order no.: 6161 000

Accessory: Shelf TH 15

Technical specifications

Temperature sensor: PT 100

Air circulation: 240 m³/h

Electrical supply: 230 V , 50 / 60 Hz (115 V version on request)

Dimensions (w x d x h): 510 x 450 x 455 mm

Weight: 21 kg