

Miniature Shaker KM 2 Akku



Description

Well suited for operation in refrigerators or incubators.

Basic Equipment

Basic device incl. basic platform and rubber mat,
with battery charger / power unit

Loading Capacity

Erlenmeyer flasks 10 ml 20 pieces
Erlenmeyer flasks 50 ml 12 pieces

Details

Motion: orbital
Shaking platform: 200 x 295 mm
Max. load: 2 kg
Shaking speed: 30 - 420 rpm, steplessly adjustable
Stroke: 8 mm
Runtime Battery: min. 24h
Battery charging time: 5h
Order No.: 6135 000

Accessory: Universal tray KM Mini
Universal tray KM
Clamping strip KM

Technical specifications

Electrical supply battery charger: 100 - 240 V, 50 / 60 Hz
Internal electrical supply: 2 batteries 24 V / 2,2 Ah
Enclosure protection: IP 21
Heat emission: approx. 4 W
Ambient temperature: 5°C to 40°C
Relative humidity: ~ 85 %
Dimensions (w x d x h): 220 x 410 x 130 mm
Weight: 9,5 kg