

Three-Dimensional Shaker TL 10



Description

Slow, three-dimensional motion

Basic Equipment

Basic device incl. basic tray and rubber mat

Details

Motion:	three-dimensional orbital motion
Shaking platform:	340 x 340 mm
Max. load:	10 kg
Shaking speed:	1 - 23 rpm, steplessly adjustable
Angle of tilt:	0° to 9° einstellbar
Timer:	0 - 120 min / continuous
Order No.:	6117 000
Accessory:	Multi-storey rack system WS / TL

Technical specifications

Electrical supply:	230 V or 115 V, 50 / 60 Hz please indicate in case of order
Enclosure protection:	IP 21

Heat emission: approx 7 - 10 W
Ambient temperature:: 5°C to 50°C
Relative humidity: ~ 85 %
Dimensions (w x d x h): 355 x 450 x 205 mm
Weight: 18,5 kg