

Universal Shaker SM 30 C control



Description

Programmable shaking tasks for the mixing of difficult media

Basic Equipment

Basic device with shaking plate and rubber mat, **without rack system**

With digital display and programmable control

Loading Capacity

Erlenmeyer flasks 100 ml	42 pieces
Erlenmeyer flasks 250 / 300 ml	20 pieces
Separating funnels 250 ml	6 pieces
Separating funnels 1000 ml	4 pieces
Test tube racks	5 pieces

Option:

[Control interface](#)

[USB interface](#)

Details

Motion: can be switched from orbital to reciprocating

Shaking platform: 560 x 400 mm

Max. load: 30 kg

Shaking speed: 15 - 300 rpm, steplessly adjustable
Stroke: 26 mm
Runtime: programmable / continuous
Order No: 6104 000
Accessory: Rack System Combifix SM A
Rack System Combifix SM B
Rack System Combifix SM C
Universal tray SM
2-storey top frame SM
Incubator Hood TH 30

Technical specifications

Electrical supply: 230 V or 115 V, 50 / 60 Hz please indicate in case of order
Enclosure protection: IP 21
Heat emission: approx. 20 - 30 W
Ambient temperature: 5°C to 50°C
Relative humidity: ~ 85 %
Dimensions (w x d x h): 680 x 610 x 160 mm
Weight: 33 kg